


Wound healing requires a holistic approach that addresses wound bed preparation as well as treating the underlying cause. Successful local wound management requires:

- Removal of devitalized tissue
- Management of bioburden
- Maintenance of moisture balance



The Wound Bed Preparation Paradigm³



Hydrofera Blue is your choice to overcome key barriers to healing


Hydrofera Blue dressings can help address all three components. Helping you create an environment where wounds can heal.¹⁻³

Overcoming Key Barriers to Healing



Choose Hydrofera Blue CLASSIC	
Polyvinyl Alcohol Foam	
Wound Characteristics	<ul style="list-style-type: none">• Deep, tunneling or flat• With or without devitalized tissue / slough & debris• Moderate to heavy exudate
Clinical Objectives	<ul style="list-style-type: none">• Debride• Manage bioburden• Help maintain moisture balance• Protect (cover) or pack

Choose Hydrofera Blue READY	
Polyurethane Foam	
Wound Characteristics	<ul style="list-style-type: none">• Flat• Moderate to heavy exudate
Clinical Objectives	<ul style="list-style-type: none">• Manage bioburden• Help maintain moisture balance• Up to 7 day wear time



Indications for use: Pressure ulcers, Diabetic ulcers, Venous stasis ulcers, Arterial ulcers, Superficial burns, Donor sites, Skin grafts, Post surgical incisions, Traumatic wounds, Skin tears, Abrasions, Lacerations

Choose Blue...
Safe, Natural Debridement

Devitalized tissue can slow wound healing and may increase the risk of infection.

- Autolytic debridement is known to be a safe and selective form of debriding wounds.
- Multiple authors have described autolytic debridement with use of Hydrofera Blue CLASSIC dressings.
- In a porcine eschar study, a statistically significant reduction in eschar was seen with use of Hydrofera Blue CLASSIC dressing.



Safe, Effective, Bioburden Management

Managing bioburden helps the body's immune system work to restore natural bacterial balance. Well-known antibacterials methylene blue and gentian violet work to help bind bacteria and harmful endotoxins to the dressings.

- Provides broad-spectrum antibacterial activity and acts as protective cover for wounds at risk of infection.
- Is non-cytotoxic.
- Does not inhibit growth factor activity.


Effective against the broad-spectrum of microorganisms commonly found in wounds.					
Bacteria					Yeast
<ul style="list-style-type: none">Methicillin-resistant Staphylococcus aureus (MRSA)Staphylococcus aureusStaphylococcus epidermidis (coagulase negative)	<ul style="list-style-type: none">Pseudomonas aeruginosaPseudomonas fluorescensEscherichia coliEnterococcus faecalis	<ul style="list-style-type: none">Vancomycin-resistant Enterococcus (VRE)Streptococcus pyogenesKlebsiella pneumoniaeProteus mirabilis	<ul style="list-style-type: none">Proteus vulgarisEnterobacter aerogenesBacillus subtilisYersinia enterocolitica	<ul style="list-style-type: none">Serratia marcescensStreptococcus pneumoniae*Enterobacter cloacae*	<ul style="list-style-type: none">Candida albicansCandida kruseiCandida glabrata

* Tested on Hydrofera Blue Ready dressing only

Optimal Moisture Balance

Inadequate moisture balance may be detrimental to wound healing. Hydrofera Blue dressing's micro-pore structure is specifically designed for managing exudate.

- Rapidly wicks moisture away from the wound.
- Helps maintain moisture balance at the wound surface.
- Is highly absorptive, well suited for moderate to heavy exuding wounds.



Ordering Information

Hydrofera Blue Antibacterial Foam Dressings

Product No	Product Size	Qty/Box
Hydrofera Blue CLASSIC Foam Dressings		
HB2214	5 cm x 5 cm (2" x 2")	10
HB4414	10.2 cm x 10.2 cm (4" x 4")	10
HB6614	15.2 cm x 15.2 cm (6" x 6")	10
HBT0906	9 mm Tunneling Dressing (1.2 g)	10 pair
HBRF2650	6.4 cm (2.5") diameter Ostomy Dressing	10

Hydrofera Blue CLASSIC Heavy Drainage Foam Dressings		
HBHD4450	10.2 cm x 10.2 cm x 1.27 cm (4" x 4" x 0.5")	5
HBHD6675	15.2 cm x 15.2 cm x 1.9 cm (6" x 6" x 0.75")	5

NEW Hydrofera Blue READY Foam Dressings		
HBRS2520	6.4 cm x 6.4 cm (2.5" x 2.5")	10
HBRS4520	10.2 cm x 12.7 cm (4" x 5")	10
HBRS8820	20.3 cm x 20.3 cm (8" x 8")	10

- Applewhite AJ, Attar P, Linden B, Stevenson Q. Gentian violet and methylene blue polyvinyl alcohol foam antibacterial dressing as a viable form of autolytic debridement in the wound bed. Surg Technol Int. 2015 May;26:65-70.
- Coutts PM, Ryan J, Sibbald RG. Case series of lower-extremity chronic wounds managed with an antibacterial foam dressing bound with gentian violet and methylene blue. Adv Skin Wound Care. 2014; 27(Suppl 1):9-13.
- Sibbald RG, Elliott JA, Ayello EA, Somayaji R. Optimizing the Moisture Management Tightrope with Wound Bed Preparation 2015. Adv Skin Wound Care. 2015 Oct;28(10):466-76
- Strohal, R., Apelqvist, J., Dissemmond, J. et al. EWMA Document: Debridement. J Wound Care. 2013; 22 (Suppl. 1): S1-S52.
- Edwards K. New twist on an old favorite: gentian violet and methylene blue antibacterial foam dressings. Adv Wound Care (New Rochelle). 2015. In Press.
- Ramsay S, Jelf C, Aust D, Roche E. Comparison of the effects of antimicrobial foam dressings on cell viability, proliferation, and growth factor activity. Poster presented at: World Union of Wound Healing Societies Congress; June 4-8, 2008, Toronto, Canada.
- Internal data on file.

Manufactured by
Hollister Incorporated
2000 Hollister Drive
Libertyville, IL 60048 USA

Distributed in Canada by
Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3

www.hollister.com

For product questions and samples please call Customer Service
at 1-800-263-7400 or your Hollister representative

The Hollister logo is a trademark of Hollister Incorporated.
Hydrofera Blue is a trademark of Hydrofera, LLC.
All other trademarks are the property of their respective owners.
©2015 Hollister Limited
922783-CA



Choose Blue.

Hydrofera Blue dressings

New Expanded Portfolio



Overcome key barriers to healing